

T-CUP URINARY TEST

Urine Screening Drug Test

- ◆ Approved by Health Canada
- ◆ Approved by the FDA (Food and Drug Administration)
- ◆ Easy to use and interpret
- ◆ Results in only 5 minutes
- ◆ Temperature indicator of the sample
- ◆ Secure seal




Detects 14 drugs and drug families

1. Amphetamines (AMP)
2. Tricyclic antidepressants (TCA)
3. Barbiturates (BAR)
4. Benzodiazepines (BZO)
5. Buprenorphine (BUP)
6. Cannabis (THC)
7. Cocaine (COC)
8. Ecstasy (MDMA)
9. Methadone (MTD)
10. Metamphetamine (mAMP)
11. Opiates (OPI)
12. Oxycodone (OXY)
13. Phencyclidine (PCP)
14. Propoxyphene (PPX)

Detects 3 types of adulterants


1. pH modifiers (pH)
2. Oxidizing and bleaching agents (Creatinine, CR)
3. Density modifiers (SG)



T-CUP URINARY TEST

Urine Screening Drug Test


INSTRUCTIONS



IMPORTANT: Refer to the front side to view a presumptive positive result for a substance (only 1 line in the control area) and no line in the test area.



1
The person to be tested must fill the container with a minimum of 40 ml of urine. Refer to the minimum filling level.



2
Close the lid and place the container on a flat surface and watch the temperature gauge. *For comparison purposes, it is advisable to take the body temperature using a thermometer.*



3
Wait 5 minutes. Remove the front label and read the result immediately.

Refer to the adulteration color chart for: creatinine, CR; Specific gravity, SG; pH, PH to validate that the sample displays normal values.

Thresholds and times of detection of substances in the urine with the use of the T-Cup

Substance	Abbreviation	Threshold	Minimum Detection time	Maximum Detection time
Amphetamines	AMP	1000 ng/ml	2 to 7 hours	2 to 4 days
Tricyclic Antidépresseurs	TCA	1000 ng/ml	8 to 12 hours	2 to 7 days
Barbiturates	BAR	300 ng/ml	2 to 4 hours	1 to 3 weeks
Benzodiazépines	BZO	300 ng/ml	2 to 7 hours	1 to 4 days
Buprenorphine	BUP	10 ng/ml	2 to 7 hours	2 to 3 days
Cannabis	THC	50 ng/ml	2 hours	Up to 40 days
Cocaine	COC	150 ng/ml	1 to 74 hours	2 to 4 days
Ectasy	MDMA	500 ng/ml	2 to 7 hours	2 to 4 days
Methadone	MTD	300 ng/ml	3 to 8 hours	1 to 3 days
Methamphetamine	mAMP	1000 ng/ml	2 to 7 hours	2 to 4 days
Opiates	OPI	2000 ng/ml	2 hours	2 to 3 days
Oxycodone	OXY	100 ng/ml	1 to 3 hours	1 to 2 days
Phencyclidine	PCP	25 ng/ml	4 to 6 hours	7 to 14 days
Propoxyphène	PPX	300 ng/ml	2 to 7 hours	2 to 3 days

CAUTION: Please read this document carefully before administering this test. The test can't be used to determine the amount of drug consumed, the level of intoxication, or determine the frequency of consumption.

TO ORDER:
alco prevention.com ♦ 1-888-863-8660

